

## Biology Lab. 2 - Enzyme assay -

	Student ID	Name
Experimenter		
Exp. partner		
Exp. period		

### Result 1: Protein concentration in blood plasma

Determine the protein concentration in provided blood plasma by measuring the absorbance

Table 1: Absorbance of blood plasma

Wave length/nm	280
Absorbance	

This value is obtained from \_\_\_\_\_ times diluted sample of blood plasma.

Protein concentration can be expressed by the following equation;

$$C = K \times A$$

C: protein concentration in blood plasma

K: conversion factor

A; Absorbance of sample

[note] The conversion factor  $K$  is 1.15.

Thus, protein concentration in undiluted blood plasma is

\_\_\_\_\_ mg/mL.

**Result 2: Measurement of enzyme activity**

**Table 2: Measurement of enzyme activities using undiluted substrate and undiluted blood plasma**

Time/min	Absorbance ( $A_{340}$ )
0.5	
1	
1.5	
2	
2.5	
3	
3.5	
4	
4.5	
5	

**Result 3: The effect of substrate concentration on enzyme activity**

Based on the results of Table 2, blood plasma is diluted \_\_\_\_\_ times not to make the difference in absorbance within 1 min exceed 0.050.

**Table 3: Changes in absorbance ( $A_{340}$ ) at different concentrations of the substrate with constant blood plasma concentration**

Time/min	Substrate concentration							
	Undiluted 1 mM	Conc-1	Conc-2	Conc-3	Conc-4	Conc-5	Conc-6	Conc-7
0.5								
1								
1.5								
2								
2.5								
3								
3.5								
4								
4.5								
5								

Based on the results of Table 3, make a graph with incubation time on transverse axis and absorbance ( $A_{340}$ ) on vertical axis (Figure 1).

[note] Add trend line and show the equation to the graph.

Calculate the reaction rate ( $v$ ) from absolute value of the slope of the curve (change in absorption per min) in each substrate concentration (Table 4).

The molar absorptivity of NADH at 340 nm is  $6.22 \times 10^3 \text{ M}^{-1}\text{cm}^{-1}$  for a cuvette with an optical path length of 1 cm.

[note] M = mole/L

**Table 4: Changes in absorbance and reaction rate in different substrate concentrations**

Substrate conc. ([S]) (mM)	1							
Change in absorption per min ( $\text{min}^{-1}$ )								
Reaction rate ( $v$ ) ( $\mu\text{M}/\text{min}$ )								

On the basis of these results in Table 4, make a graph (Figure 2) with substrate concentration ([S]) on transverse axis and the reaction rate ( $v$ ) on vertical axis. Also make a Lineweaver-Burk plot (Figure 3) by taking double-reciprocal values ( $1/[S]$  as a transverse axis and  $1/v$  as vertical axis).

[note] Substrate concentration can be calculated with undiluted substrate concentration and the added volume to the reaction.

(example) two times dilution

$$30 \text{ mM} \div 2 \times (0.1 \text{ mL} \div 3 \text{ mL}) = 0.5 \text{ mM}$$

[note] Add the trend line to Figure 2. Draw the trend line to the x-intercept and add the equation to Figure 3.

Based on the equation in Figure 3, the values of  $K_m$  and  $V_{max}$  of enzyme in blood plasma can be calculated as follows;

$$K_m = \underline{\hspace{2cm}} \text{ mM}$$

$$V_{max} = \underline{\hspace{2cm}} \mu\text{M}/\text{min}$$

Details of the calculation (You may add extra pages if needed.)

## Result 4: The effect of blood plasma concentration on enzyme activity

**Table 5: Changes in absorbance ( $A_{340}$ ) at different concentrations of the blood plasma with constant substrate concentration**

Time/ min	Un-diluted (mg/mL)	Dilution rate: 1/2	Dilution rate:	Dilution rate:	Dilution rate:	Dilution rate:	Dilution rate:
0.5							
1							
1.5							
2							
2.5							
3							
3.5							
4							
4.5							
5							

Based on Table 5, make a graph (Figure 4) with time on transverse axis and absorbance ( $A_{340}$ ) on vertical axis.

[note] Add trend line and corresponding equation to Figure 4.

Calculate the reaction rate ( $v$ ) from absolute value from the slopes of the equations in each substrate concentration (Table 6).

[note] Convert dilution rates into amounts of blood plasma (E).

**Table 6: Changes in absorbance and reaction rate in different blood plasma concentrations**

Dilution rate	1	1/2					
Amount of blood plasma (E) (mg)							
Change in $A_{340}$ ( $\text{min}^{-1}$ )							
Reaction rate ( $v$ ) ( $\mu\text{M}/\text{min}$ )							

Based on Table 6, make a graph (Figure 5) with the amount of blood plasma (E) on transverse axis and the reaction rate (v) on vertical axis.

[note] The amount of blood plasma can be calculated with undiluted blood plasma concentration (from Result 1) and the added volume to the reaction.

(example) If you diluted 50 mg/mL two times,  
 $50 \text{ mg/mL} \div 2 \times 0.1 \text{ mL} = 2.5 \text{ mg}$

[note] Add the trend line and the equation to Figure 5.